




Product Specification
SafeSeal reaction tube, 2 ml, PP, PCR
Performance Tested, Low protein-binding

Page 1



Product description

Order number	72.695.600
Product description	SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low protein-binding, 50 piece(s)/minigrip bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
SafeSeal reaction tube, 2 ml, PP, PCR
Performance Tested, Low protein-binding

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	800
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	200
Piece (s) / outer case	800
Piece(s) / pallet	62400
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0011 kg
Weight of case	1.56 kg
EAN of inner box	4038917379281
EAN case	4038917379274

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT