



Product description

Order number	72.704.700
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low DNA-binding, 50 piece(s)/minigrip bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	0,5 ml
Centrifugation max (RCF)	30000 x g

Size

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1200
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	300
Piece (s) / outer case	1200
Piece(s) / pallet	93600
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0005 kg
Weight of case	1.28 kg
EAN of inner box	4038917350037
EAN case	4038917350020

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

