



Product description

| | |
|---------------------|--|
| Order number | 72.706.600 |
| Product description | SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low protein-binding, 50 piece(s)/minigrip bag |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Cap | SafeSeal cap |
| Volume of work | 1,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 39 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
SafeSeal reaction tube, 1.5 ml, PP, PCR
Performance Tested, Low protein-binding

Page 2

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 800 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 200 |
| Piece (s) / outer case | 800 |
| Piece(s) / pallet | 62400 |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.001 kg |
| Weight of case | 1.51 kg |
| EAN of inner box | 4038917264174 |
| EAN case | 4038917264167 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT